

Stage 3 Upgrades

Improvements over the State 2 Control

- Touch pad control versus dip switches allows simplified user entry of parameters
- Express cycle available on Vacumatic units
- RTD sensor versus thermistor for the vessel temperature and pressure control
- Automated utility shutdown
- Duplicate print capability to assist with record keeping
- Simplified piping configuration
- Simplified wiring using Control Interface Connector Board allows greater reliability and easier troubleshooting
- Water in jacket sensor
- Touchpad calibration of unit versus manual potentiometers
- Power door diagnostics available in service mode